

Activity Assay for Rat PAI-1

For Research Use Only

INTRODUCTION

Plasminogen Activator Inhibitor-1 (PAI-1) is a glycoprotein and member of the serine proteinase inhibitor (serpin) superfamily. PAI-1 is the primary inhibitor of tissue-type plasminogen activator (tPA) and the urokinase-type plasminogen activator (uPA). This inhibition exhibits antiproteolytic properties that can lead to myocardial infarction and thromboembolic disease with elevated levels of PAI-1. Additionally, PAI-1 is thought to play a role in the function of tissue remodeling and tumor metastasis.

PRINCIPLES OF PROCEDURE

This is an ELISA (Enzyme-Linked Immunosorbent Assay) for the quantitative analysis of PAI-1 levels in biological fluid. This test kit operates on the basis of sandwich ELISA where active PAI-1 complexes with uPA and is quantified with the use of an HRP labeled secondary antibody.

The functional or active PAI-1 binds to the uPA coated on the well of the microtiter plate. The latent and complexed forms of PAI-1 will not bind and are discarded at a later washing step. Next, PAI-1 primary antibody is added to the wells binding to the captured PAI-1 on the microtiter plate. HRP conjugated secondary antibody is then added for detection of the active PAI-1. Quantitative test results are obtained by the measure and comparison of the sample and standard absorbance readings.

MATERIALS PROVIDED

Component	Contents	Quantity	Storage	Cat. No.
Coated Plate	96-well microplate coated with Anti-PAI-1	1 plate	4°C	PI92a
Standard	Rat PAI-1 activity standard (50 ng/mL)	1 vial	4°C	PI92b
Wash Buffer	10x solution for washing plate	50 mL	4°C	PI92c
Substrate	TMB Substrate	10 mL	4°C	PI92d
Primary Antibody	Anti-rat PAI-1 polyclonal antibody (lyophilized)	1 vial	4°C	PI92e
Secondary Antibody	Anti-rabbit HRP conjugated antibody	1 vial	4°C	PI92f

MATERIALS NEEDED BUT NOT PROVIDED

- 1 N H₂SO₄
- TBS Buffer (see Reagent Preparation)
- Blocking Buffer (see Reagent Preparation)
- DI Water
- Microplate reader with 450 nm filter
- Microplate shaker with uniform horizontal circular movement up to 300 rpm
- Precision pipettes that range from 10 µL-1000 µL and disposable tips

STORAGE

1. Store the kit and all of its components at 4°C before use.

2. If not using the entire plate at once, prepare only the appropriate amount of Primary Antibody and PAI-1 Standard. The remaining stock solutions should be frozen and stored at -70°C . Primary Antibody should be used within two weeks. All other components should remain refrigerated.
3. Store unused portions of the microplate in a pouch with a desiccant at 4°C .

PROCEDURAL NOTES

1. This assay should be run at room temperature.
2. Use aseptic technique when opening and dispensing reagents.
3. This kit is designed to work properly as provided and instructed. Additions, deletions or substitutions to the procedure or reagents are not recommended, as they may be detrimental to the assay.
4. To minimize error due to handling, wipe the exterior bottom of the microplate wells with a lint-free paper towel.

SAMPLE COLLECTION, STORAGE, AND PREPARATION

Blood should be collected in a 9:1 ratio with trisodium citrate or acidified citrate. Immediately after collection, centrifuge the samples at $3000 \times g$ for 15 minutes. The plasma must be platelet free as they can release PAI-1. Store the plasma samples on ice prior to analysis. The plasma is stable at -20°C for one month with 3 freeze-thaw cycles.

REAGENT PREPARATION

1. **TBS Buffer:** 0.10 M TRIS, 0.15 M NaCl, pH 7.4
2. **Blocking Buffer:** 3% BSA in TBS Buffer.
3. **10x Wash Buffer:** Dilute the 50 mL of concentrate to 1x with 450 mL of DI water prior to use.
4. **Standard:** Reconstitute with 3% BSA Blocking Buffer as directed on the vial and vortex gently to mix. Prepare according to include Standard Dilution Table immediately prior to use.
5. **Primary Antibody:** Reconstitute with 3% BSA Blocking Buffer as directed on the vial and vortex gently to mix. Prepare immediately prior to use.
6. **Secondary Antibody:** Dilute with 3% BSA Blocking Buffer as directed on the vial and vortex gently to mix. Prepare immediately prior to use.

ASSAY PROCEDURE


1. Add 100 μL of standards or unknowns to each well. See **Scheme I** for a suggested plate layout. Shake the plate at 300 rpm on a plate shaker for 30 minutes.
2. Wash wells according to the following wash procedure:
 - a. Remove contents of the plate by inversion into an appropriate disposal device.
 - b. Tap out remaining contents of the plate onto a lint free paper towel.
 - c. Add 300 μL of Wash Buffer to each well.
 - d. Let stand for 2-3 minutes.
 - e. Remove contents of the plate by inversion into an appropriate disposal device.
 - f. Repeat procedure 2 more times and proceed to step "g".
 - g. Tap out the remaining contents of the plate onto a lint free paper towel and proceed to step 6.
3. Add 100 μL of diluted Primary Antibody to each well. Shake plate at 300 rpm for 30 minutes.
4. Wash wells according to step 2.
5. Add 100 μL of the diluted Secondary Antibody to each well. Shake plate at 300 rpm for 30 minutes.
6. Wash wells according to step 2.
7. Add 100 μL of TMB Substrate to each well and incubate for 2-10 minutes with shaking.
8. Stop the reaction with 50 μL per well of 1 N H_2SO_4 and read plate at 450 nm.

Scheme I:

	1	2	3	4	5	6	7	8	9	10	11	12
A	S ₀	S ₁	S ₂	S ₃	S ₄	S ₅	S ₆	S ₇	S ₈	S ₉	S ₁₀	U ₁
B	S ₀	S ₁	S ₂	S ₃	S ₄	S ₅	S ₆	S ₇	S ₈	S ₉	S ₁₀	U ₁
C	U ₂	U ₃	U ₄	U ₅	U ₆	U ₇	U ₈	U ₉	U ₁₀	U ₁₁	U ₁₂	U ₁₃
D	U ₂	U ₃	U ₄	U ₅	U ₆	U ₇	U ₈	U ₉	U ₁₀	U ₁₁	U ₁₂	U ₁₃
E	U ₁₄	U ₁₅	U ₁₆	U ₁₇	U ₁₈	U ₁₉	U ₂₀	U ₂₁	U ₂₂	U ₂₃	U ₂₄	U ₂₅
F	U ₁₄	U ₁₅	U ₁₆	U ₁₇	U ₁₈	U ₁₉	U ₂₀	U ₂₁	U ₂₂	U ₂₃	U ₂₄	U ₂₅
G	U ₂₆	U ₂₇	U ₂₈	U ₂₉	U ₃₀	U ₃₁	U ₃₂	U ₃₃	U ₃₄	U ₃₅	U ₃₆	U ₃₇
H	U ₂₆	U ₂₇	U ₂₈	U ₂₉	U ₃₀	U ₃₁	U ₃₂	U ₃₃	U ₃₄	U ₃₅	U ₃₆	U ₃₇

CALCULATIONS

1. Average the O.D. values for each pair of duplicate wells.
2. Plot a standard curve using the average O.D. versus the standard concentration.
3. Fit a straight line through the points using a linear fit procedure.
4. Determine the concentration of each unknown using the equation from the standard curve.

Typical Standard Curve:**PERFORMANCE CHARACTERISTICS**

Sensitivity: 0.011 ng/mL

REFERENCES

1. Thorton, A., *et al.*, (1994) *Cancer Research* 54(5):1337-1343
2. Hamdan, A., *et al.*, (1996) *Circulation* 93:1073-1078
3. DeYoung, M., *et al.*, (2001) *Circulation* 104:1971-1972
4. Hasenstab, D., *et al.*, (2000) *Arterioscler Thromb Vasc Biol* 20:853-859
5. Hagiwara, H., *et al.*, (2003) *Thromb Res* 111(4-5):301-309
6. Barnes, J., *et al.*, (1995) *J. Histochem Cytochem* 43(9):895-905
7. van Leeuwen, R., *et al.*, (1994) *Circulation* 90:362-368
8. Thu-Hoa, N., *et al.*, (1998) *Thromb Haemost* 79:808-812

DISCLAIMER

This information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Oxford Biomedical Research, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See catalog for additional terms and conditions of sale.

ORDERING INFORMATION

For additional kits or a complete catalog please call 800-692-4633.

TECHNICAL SUPPORT

If you need technical information or assistance with assay procedures, call our Technical Support Department at 800-692-4633 or 248-852-8815. Our staff will be happy to answer your questions about this or any other product in the Oxford Biomedical line.

GUARANTEE AND LIMITATION OF REMEDY

Oxford Biomedical Research, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the material in this ELISA kit, except that these materials and this kit will meet our specifications at the time of delivery. Buyer's remedy and Oxford Biomedical Research, Inc.'s sole liability hereunder is limited to, at Oxford Biomedical Research, Inc.'s option, refund of the purchase price of, or the replacement of, material that does not meet our specification. By acceptance of our products, Buyer indemnifies and holds Oxford Biomedical Research, Inc. harmless against, assumes all liability for the consequence of its use or misuse by the Buyer, its employees, or others. Said refund or replacement is conditioned of Buyer notifying Oxford Biomedical Research, Inc. within (30) days of the receipt of product. Failure of Buyer to give said notice within said thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

Oxford Biomedical Research, Inc.
P.O. Box 522
Oxford, MI 48371 U.S.A.

Orders: 800-692-4633
Technical Service: 248-852-8815
Fax: 248-852-4466
E-mail: info@oxfordbiomed.com

Made in the U.S.A.