

Microplate Assay for Lipid Peroxidation


For Research Use Only

INTRODUCTION

Lipid peroxidation is a well-established mechanism of cellular injury in both plants and animals, and is used as an indicator of oxidative stress in cells and tissues. Lipid peroxides are unstable and decompose to form a complex series of compounds including reactive carbonyl compounds. Polyunsaturated fatty acid peroxides generate malondialdehyde (MDA) and 4-hydroxyalkenals (HAE) upon decomposition. The measurement of MDA has been used as an indicator of lipid peroxidation (1).

PRINCIPLES OF PROCEDURE

This assay is based on the reaction of a chromogenic reagent, N-methyl-2-phenylindole (R1), with MDA at 45°C. One molecule of MDA reacts with 2 molecules of Reagent R1 to yield a stable chromophore with maximal absorbance at 586 nm.



MDA : R = OH

4-hydroxyalkenal : R= hydroxyalkyl

Max = 586 nm

MATERIALS PROVIDED

Component	Description	Volume	Storage
Reagent R1	N-methyl-2-phenylindole in Acetonitrile	13 mL	4°C
MDA Standard	1,1,3,3-Tetramethoxypropane (TMOP) in Tris-HCl	200 µL	4°C
Diluent	Ferric Iron in Methanol	5.5 mL	4°C
Microplate	96-well Microplate	1 plate	4°C

MATERIALS NEEDED BUT NOT PROVIDED

1. Spectrophotometer plate reader with a 586 nm filter (580 nm or 590 nm filters are also acceptable)
2. Water bath at 45°C
3. 37% HCl (optional – used for measuring MDA only)
4. Butylated hydroxytoluene (BHT) (used in sample collection)
5. Acetonitrile
6. Microcentrifuge
7. Polypropylene microcentrifuge tubes
8. Deionized Water (dH₂O)
9. Adjustable micropipettes (10 – 1000 µL) and tips

Optional:

1. Methansulfonic Acid (MSA)

STORAGE

1. The reagents are stable until the indicated kit expiration date if handled and stored properly.
2. When not in use, place the bottles at 4°C.

WARNINGS AND PRECAUTIONS

1. Use aseptic technique when opening and dispensing reagents.
2. In case of accidental exposure of skin, mucous membranes or eyes to R1 or R2 reagents, thoroughly wash the exposed area with water.
3. This kit is designed to work properly as provided and instructed. Additions, deletions or substitutions to the procedure or reagents are not recommended, as they may be detrimental to the assay.

PROCEDURAL NOTES

1. Do not leave the reagent bottles open. Replace the caps as soon as the desired volume is removed.
2. Do not allow the capped reagent bottles to sit at room temperature for long periods of time.
3. 4-hydroxyalkenals (HAE) plus MDA can be measured by substituting methanesulfonic acid for HCl in **ASSAY PROCEDURE** Step 3 below.
4. To minimize error due to handling, wipe the exterior bottom of the microplate wells with a lint-free paper towel.

REAGENT PREPARATION

1. **100x BHT Stock Solution:** 0.5 M in Acetonitrile.
2. **Reagent R1:** Dilute Reagent R1 3:1 with Diluent (i.e. 12 mL R1 + 4 mL Diluent). Prepare immediately before use.
3. **37% HCl:** 12 N HCl – Do not dilute prior to use in the assay.
4. **Sample Blank:** 75% Acetonitrile/ 25% Diluent. Add 650 μ L to a microcentrifuge tube. Skip assay steps 1 and 2. The acid addition and the incubation steps are carried out on this blank.
5. **20 μ M MDA Standard Stock:** Dilute the 10 mM MDA Standard 1:500 in dH₂O (i.e. 20 μ L 10 mM MDA + 9.98 mL dH₂O). Prepare immediately before use.

SAMPLE PREPARATION

Sample Stability

Unless assayed immediately, samples should be frozen at -70°C to prevent loss of MDA (3,4) and prevent new sample oxidation. Samples should not be stored at -20°C. Once thawed from -70°C storage for assay, the sample should not be refrozen.

Samples should be protected from light to avoid photooxidation.

Oxidation Prevention

We recommend adding BHT at a final concentration of 5 mM in the buffer prior to homogenization of tissue or cells. If no antioxidant is added, new lipid peroxidation can occur during homogenization and biased values will result (2).

Plasma or Serum

The amount of free MDA in normal plasma or serum is at or below the limit of detection of this assay.

Cell Culture (8,9)

Cells cultured in serum containing medium should be washed several times to remove serum components prior to homogenization. Since MDA exists as the water-soluble enolate anion at physiological pH, much of the MDA generated from lipid peroxidation in cell culture may be in the culture medium.

1. Remove cells using a rubber policeman. Lysis buffers have a high potential of interfering in this assay. Cells should be washed well in ice-cold 20 mM PBS or Tris buffer, pH 7.4, and resuspended in the same buffer. Researchers should determine the optimal number of cells to use in this assay, but a recommended starting point is 5×10^7 cells per ml.
2. Lyse cells by sonication, homogenization, or freeze-thaw cycles. To prevent sample oxidation during preparation, lysis should be done in the presence of 10 μ L 0.5 M BHT per 1 ml of cell homogenate.
3. After homogenization, follow steps 5-7 of the tissue homogenization procedure.

Tissue Homogenates (6,7)

Sample homogenates should be made as concentrated as possible. The concentration of protein in the homogenate should be determined. It is recommended that 0.2 mL of a homogenate containing 15-60 mg/mL protein be assayed for initial studies in a previously untested biological sample.

1. If necessary, remove blood *in situ* by perfusion or *in vitro* by rinsing with ice-cold isotonic saline (0.9% NaCl).
2. Weigh tissue. A reasonable amount to start with is 1 g tissue per 10 mL of buffer.
3. Prepare tissue homogenate in ice-cold 20 mM PBS or Tris buffer, pH 7.4. Other buffers may be used, but the researcher should confirm non-interference in the assay by measuring the MDA and/or 4-HNE standards diluted in the chosen buffer.
4. Add 10 μ L 0.5 M BHT Stock Solution per 1 mL of tissue homogenate to prevent sample oxidation. A precipitate is expected and will not affect the outcome of the assay as it is removed by centrifugation.
5. After homogenization, centrifuge at 3,000 x g and 4°C for 10 minutes to remove large particles.
6. Remove an aliquot of the sample for protein determination.
7. Freeze the sample immediately at -70°C or keep on ice prior to testing. Test 0.2 mL of the homogenate in the assay.

STANDARD CURVE PREPARATION

Malondialdehyde is provided as an acetal because the aldehyde itself is not stable. The acetal (TMOP) is hydrolyzed during the acid incubation step at 45°C, which will generate MDA.

Please see the **Reagent Preparation** section for preparing the 20 μ M MDA Standard Stock.

Table 1: Standard Curve Preparation

Standard	MDA Conc. (μ M)	Vol. of dH ₂ O (μ L)	Vol. of 20 μ M MDA Stock (μ L)
S ₀	0	200	-
S ₁	0.5	195	5
S ₂	1.0	190	10
S ₃	2.5	175	25
S ₄	5.0	150	50
S ₅	10.0	100	100
S ₆	15.0	50	150
S ₇	20.0	-	200

ASSAY PROCEDURE

1. Add 140 μL of Standards or Samples to a microcentrifuge tube.
2. Add 455 μL of diluted Reagent R1 to each tube and vortex.
3. Add 105 μL 37% HCl (12 N HCl) to each tube and mix well.
4. Incubate at 45°C for 60 minutes.
5. Centrifuge samples at 15,000 x *g* for 10 minutes to obtain a clear supernatant.
6. Transfer the 3 x 150 μL of the supernatant to the microplate and read at 586 nm. See Scheme I for a sample plate layout. The color is stable for at least an hour at room temperature, or 2 hours at 4°C when stored in the dark.

Scheme I: Sample Plate Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	S ₀	S ₀	S ₀	BLK	BLK	BLK	U ₈	U ₈	U ₈	U ₁₆	U ₁₆	U ₁₆
B	S ₁	S ₁	S ₁	U ₁	U ₁	U ₁	U ₉	U ₉	U ₉	U ₁₇	U ₁₇	U ₁₇
C	S ₂	S ₂	S ₂	U ₂	U ₂	U ₂	U ₁₀	U ₁₀	U ₁₀	U ₁₈	U ₁₈	U ₁₈
D	S ₃	S ₃	S ₃	U ₃	U ₃	U ₃	U ₁₁	U ₁₁	U ₁₁	U ₁₉	U ₁₉	U ₁₉
E	S ₄	S ₄	S ₄	U ₄	U ₄	U ₄	U ₁₂	U ₁₂	U ₁₂	U ₂₀	U ₂₀	U ₂₀
F	S ₅	S ₅	S ₅	U ₅	U ₅	U ₅	U ₁₃	U ₁₃	U ₁₃	U ₂₁	U ₂₁	U ₂₁
G	S ₆	S ₆	S ₆	U ₆	U ₆	U ₆	U ₁₄	U ₁₄	U ₁₄	U ₂₂	U ₂₂	U ₂₂
H	S ₇	S ₇	S ₇	U ₇	U ₇	U ₇	U ₁₅	U ₁₅	U ₁₅	U ₂₃	U ₂₃	U ₂₃

CALCULATIONS

1. Using the standard data, calculate the net A_{586} for each standard by subtracting the blank (S_0) value from each of the standard A_{586} values. Plot the net A_{586} against MDA concentration, and perform a linear regression analysis:


$$[A_{586}] = m[\text{MDA}] + b$$

2. Calculate the concentration of analyte in each unknown from the net A_{586} of the sample. If a sample blank was required, subtract the absorbance of the SB (A_{SB}) from the net sample absorbance:

$$[\text{MDA}] = \frac{A_{586} - b}{m} \cdot \text{df}$$

Where: $[\text{MDA}] = \mu\text{M}$ concentration of MDA in the sample $b = \text{intercept}$
 $A_{586} = \text{Net absorbance at 586 nm of the sample}$ $\text{df} = \text{dilution factor}$
 $m = \text{regression coefficient (slope)}$

Figure 1: Typical Standard Curve



PERFORMANCE CHARACTERISTICS

Least Detectable Concentration

Experiments on enaldehyde standards and blanks suggest that the analytical limit of detection in a purified system will be 0.1 nmol/mL final concentration (0.5 nmol/mL in the sample), corresponding to an absorbance value of approximately 0.011. The limit for biological samples will generally be higher and it is recommended that the researcher assess this parameter in their system.

Reproducibility

Experiments in which standard samples (0-20 μ M) were assayed using the same protocol over a period of 10 days established the standard error of the measurement (SEM) at less than 5%.

Limitations

1. Sucrose or fructose at concentrations of ≥ 50 mM in the sample will cause a high bias in the assay. Vitamin E, when present at ≥ 15 μ M, can cause a diminution in the values obtained for 4-HNE.
2. Although the standards in this assay will usually appear blue, the samples or blanks sometimes appear another color, such as pink or green. This is due to chromophores that form other than those producing the 586 nm peak. Ordinarily, these chromophores will not interfere with the A_{586} .
3. This assay measures only free MDA or 4-HNE in samples. The conditions of the assay do not provide for liberation of MDA bound to proteins via Schiff Base. 4-HNE is sufficiently reactive that it rapidly combines with proteins in tissues, forming stable adducts that are not liberated by heating at high temperatures in acid; as a consequence, there is very little free 4-HNE in tissue (1).
4. Normal tissues have very low levels of free MDA or 4-HNE, typically 10-100 pmol/mg protein (2,7). Assay of a 0.2 mL sample containing 10 mg of protein derived from normal tissue will give absorbance values at 586 nm of 0.01 or less in this assay. Caution must be taken not to interpret very low absorbance values as an accurate reflection of analyte concentrations in biological samples.
5. In setting up this assay for the first time on a particular biological sample, the kinetics of color development on the sample should be followed in comparison with that of the TMOP standard. The A_{586} of the sample should reach a plateau and then remain stable. If the A_{586} continues to go up after the standards have achieved a stable color, the researcher should be concerned that non-MDA reactivity (interference) is occurring in the sample.
6. In setting up this assay for the first time on a particular biological sample, a wavelength scan from 450 to 700 nm should be performed on the clarified sample reaction mixture and compared to that obtained with the TMOP standard. The lack of a peak at 586 nm or lack of reasonable definition to the sample profile compared to the standard would suggest interference in the sample.
7. If no antioxidant is added to the samples during homogenization and subsequent assay, a high bias due to new sample oxidation may result.

REFERENCES

1. Esterbauer, H., *et al.*; (1991) *Free Rad. Biol. Med.* **11**:81-128
2. Botsoglou, N.A., *et al.*; (1994) *J. Agric. Food Chem.* **42**:1931-1937
3. Carbonneau, M.A., *et al.*; (1991) *Clin. Chem.* **37**:1423-1429
4. Bull, A.W. and Marnett, L.J.; (1985) *Analyt. Biochem.* **149**:284-290
5. Liu, J., *et al.*; (1997) *Analyt. Biochem.* **245**:161-166
6. Mattson, J.P., *et al.*; (2002) *Pathophysiology* **8**:215-221
7. Calvo, J.R., *et al.*; (2001) *J. Cell Biochem.* **80**:461-70
8. Janciauskiene, S. and Bo Ahren; (2000) *Biochem. Biophys. Res. Com.* **267**:619-625
9. Montilla-Lopez, P., *et al.*; (2002) *European J. Pharmacology* **451**:237-243

DISCLAIMER

This information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Oxford Biomedical Research, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See catalog for additional terms and conditions of sale.

ORDERING INFORMATION

For additional kits or a complete catalog please call 800-692-4633.

TECHNICAL SUPPORT

If you need technical information or assistance with assay procedures, call our Technical Support Department at 800-692-4633 or 248-852-8815. Our staff will be happy to answer your questions about this or any other product in the Oxford Biomedical line.

GUARANTEE AND LIMITATION OF REMEDY

Oxford Biomedical Research, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the materials in this kit, except that these materials and this kit will meet our specifications at the time of delivery. Buyer's remedy and Oxford Biomedical Research, Inc.'s sole liability hereunder is limited to, at Oxford Biomedical Research, Inc.'s option, refund of the purchase price of, or the replacement of, material that does not meet our specification. By acceptance of our products, Buyer indemnifies and holds Oxford Biomedical Research, Inc. harmless against, assumes all liability for the consequence of its use or misuse by the Buyer, its employees, or others. Said refund or replacement is conditioned on Buyer notifying Oxford Biomedical Research, Inc. within thirty (30) days of the receipt of product. Failure of Buyer to give said notice within thirty (30) days of receipt of product shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

Oxford Biomedical Research, Inc.
P.O. Box 522
Oxford, MI 48371 U.S.A.

Orders: 800-692-4633
Technical Service: 248-852-8815
Fax: 248-852-4466
E-mail: info@oxfordbiomed.com

Made in the U.S.A.