

**Product Specification Sheet**

Citrullinated Human Fibrinogen

Product Number: NF31

Aliquot: 1 mg

Lot Number: nf31.220304

Storage: -70°C

**CONCENTRATION:** 0.75 mg/mL total protein.

**STORAGE BUFFER:** 100mM Tris pH 8.0, 140mM NaCl, 5mM Imidazole.

**STORAGE:** -70°C; **AVOID MULTIPLE FREEZE-THAW CYCLES.**

**PURITY:** >90%

**MOLECULAR WEIGHT:** Mixture.

**SOURCE:** Human Plasma.