

**Product Specification Sheet**

Description: CypExpress™ Null  
Product Number: CENull.1  
Aliquot: 1.0 Gram  
Lot Number: 160314  
Storage: -70°C

**DESCRIPTION:** CypExpress™ Null is a permeabilized and stabilized dried yeast powder preparation that has been sham infected and processed similarly to the other CypExpress™ p450 products. Its intended use is as a negative control.

**THAWING:** It is important to avoid water adsorption during thawing. Remove the CypExpress™ vial from storage and warm to room temperature prior to opening. Return unused powder to -70°C.

**STORAGE:** Store at -70°C. Although CypExpress™ preparations can be warmed to room temperature multiple times, it is best to aliquot as needed to avoid multiple freeze/thaw cycles.

**SUGGESTED PROTOCOLS:** The basic protocols attached are intended as a starting point that can be optimized by an experienced user for specific applications and samples.