

Product Specification Sheet

Anti-Factor X (Mouse) Biotin
Labeled
Product Number: CF06
Aliquot: 1.0 mg
Lot Number: 912
Storage: -70°C

DESCRIPTION:	This polyclonal sheep antibody is specific for mouse factor X and is biotin labeled.
TARGET MW:	58.9 kDa
CONCENTRATION:	5.0 mg/mL
SIZE:	1.0 mg
VOLUME:	0.114 mL
STORAGE:	-70°C – Avoid repeated freeze/thaw cycles.
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
SOURCE:	Sheep serum
FORM:	Frozen liquid
SPECIES SPECIFICITY:	Mouse
PURIFICATION:	IgG fraction purified by immobilized Protein G.

For research purposes only.