

Mouse Factor X  
Product Number: CF38  
Aliquot: 0.1 mg  
Lot Number: cf38.514  
Storage: -70°C

<b>DESCRIPTION:</b>	This mouse Factor X is isolated from fresh frozen plasma and is purified immunoaffinity chromatography.
<b>MOLECULAR WEIGHT:</b>	58.9 kDa
<b>CONCENTRATION:</b>	1.4 mg/mL
<b>SIZE:</b>	0.1 mg
<b>STORAGE:</b>	-70°C – <b>Avoid repeated freeze/thaw cycles.</b>
<b>BUFFER:</b>	5 mM MES; 150 mM NaCl; pH 6.0
<b>PURIFICATION:</b>	Affinity purified
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Mouse plasma

For Research Purposes Only.