

Product Specification Sheet

Anti-Prothrombin (Mouse) IgG
Fraction
Product Number: CT04
Aliquot: 1.0 mg
Lot Number: ct04.607
Storage: -70°C

DESCRIPTION:	This polyclonal sheep antibody is specific for mouse prothrombin/thrombin and is provided as a purified IgG fraction.
MOLECULAR WEIGHT:	72 kDa
CONCENTRATION:	5.7 mg/mL
SIZE:	1.0 mg
VOLUME:	0.175 mL
STORAGE:	-70°C – Avoid repeated freeze/thaw cycles.
IMMUNOGEN:	Purified, full length native protein from mouse.
SOURCE:	Sheep serum
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
PURIFICATION:	IgG fraction purified by immobilized Protein G
SPECIES SPECIFICITY:	Mouse, Human

For Research Purposes Only.