

## Product Specification Sheet

Rabbit Fibrinogen  
Product Number: CT17  
Aliquot: 25 mg  
Lot Number:  
Storage: -70°C

**DESCRIPTION:** This is fibrinogen purified from plasma using a multi-step precipitation process followed by SPE extraction with immunoaffinity resins.

**MOLECULAR WEIGHT:** ~330 kDa

**CONCENTRATION:** 5.6 mg/mL determined spectrophotometrically.

**STORAGE BUFFER:** 0.01 M HEPES; 0.3 M NaCl; 10 U/mL Trasylol; pH 7.4

**STORAGE:** -70°C **AVOID REPEATED FREEZE THAW CYCLES**

**PURIFICATION:** Purified using a multi-step precipitation process followed by SPE extraction with immunoaffinity resins.

**NOTES:** This product is >90% clottable and is homogenous on 4-20% SDS-PAGE gels.

**REFERENCES:**