

Product Specification Sheet

Anti-Mouse Fibrinogen Antisera
Product Number: CT19
Aliquot: 1.0 mL
Lot Number ct19.608
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as neat antiserum.
TARGET MOLECULAR WEIGHT:	340 kDa
STORAGE BUFFER:	Not Applicable
STORAGE:	-70°C
IMMUNOGEN:	Full length native mouse fibrinogen isolated from plasma using several chromatographic steps.
PURITY:	Serum
CONCENTRATION:	Whole antiserum
SPECIES SPECIFICITY:	Mouse and Rat
APPLICATIONS:	ELISA, WB, IHC-P, IHC-F