

## Product Specification Sheet

Anti-Mouse Fibrinogen HT  
Affinity Purified  
Product Number: CT27  
Aliquot: 100 µg  
Lot Number: ct27.803  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is specific for mouse fibrinogen.
<b>MOLECULAR WEIGHT:</b>	330 kDa
<b>CONCENTRATION:</b>	7.3 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Purified mouse fibrinogen from plasma.
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Affinity purified
<b>PURITY:</b>	≥ 98% assessed by visual inspection of a Coomassie® blue-stained SDS-PAGE gel.
<b>APPLICATIONS:</b>	Suitable for use in Western blot (1:10,000) and ELISA (1:100,000).
<b>SPECIES SPECIFICITY:</b>	Mouse