

**Product Specification Sheet**

Description: Mouse Thrombin  
Product Number: CT38  
Aliquot: 0.05 mg  
Lot Number: ct38.2617  
Storage: -70°C

<b>DESCRIPTION:</b>	Mouse Thrombin
<b>MOLECULAR WEIGHT:</b>	36.7 kDa
<b>CONCENTRATION:</b>	1.27 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Na Citrate; 0.2 M NaCl; 0.1% PEG-8000; pH 6.5
<b>STORAGE:</b>	-70°C
<b>SOURCE:</b>	Mouse plasma
<b>PURIFICATION:</b>	Prepared from purified prothrombin.
<b>PURITY:</b>	>95%
<b>SPECIFIC ACTIVITY:</b>	3013 NIH u/mg