

Product Specification Sheet

Description: Human Fibrinogen,
FITC Labeled
Product Number: CT53
Aliquot: 10 mg
Lot Number: ct53.typical
Storage: -70°C

DESCRIPTION:	This human fibrinogen was purified by chromatography and plasminogen has been depleted. It was then labeled with FITC at the amino terminal.
MOLECULAR WEIGHT:	330 kDa
CONCENTRATION:	10 mg/mL
STORAGE BUFFER:	20 mM Sodium Citrate; pH 7.4
STORAGE:	Store at -70°C
PURIFICATION:	Purified by multiple chromatographic steps. Plasminogen depleted.
PURITY:	>95% by SDS-PAGE analysis
USAGE:	Thaw at 37°C until completely liquid. Mix very gently before use.