

## Product Specification Sheet

Porcine Fibrinogen, FITC Labeled  
Product Number: FB04  
Aliquot: 10.0 mg  
Lot Number: fb04.616  
Storage: -70°C

<b>DESCRIPTION:</b>	Porcine fibrinogen that has been amino terminal labeled with FITC. Fibrinogen is prepared from fresh porcine plasma that has been plasminogen depleted.
<b>MOLECULAR WEIGHT:</b>	330 kDa
<b>CONCENTRATION:</b>	10.0 mg/mL
<b>STORAGE BUFFER:</b>	0.02M Sodium Citrate-HCl; pH 7.4
<b>STORAGE:</b>	Store at -70°C and protect from light.
<b>LABEL:</b>	Labeled with FITC (fluorescein isothiocyanate). 1.5 moles FITC: 1 mole protein.
<b>SPECIFIC ACTIVITY:</b>	73% Clottable.
<b>PURIFICATION:</b>	>95% pure by SDS-PAGE analysis
<b>THAWING:</b>	Fibrinogen should be thawed at 37°C without agitation. Once it is completely thawed, GENTLY mix before use. Fibrinogen should be kept at 25-37°C while using. Flash freeze aliquots.