

Product Specification Sheet

FITC Labeled Rabbit Fibrinogen
Product Number: FB05
Aliquot: 10 mg
Lot Number: fb05.411
Storage: -70°C

DESCRIPTION:	Rabbit fibrinogen amino terminal labeled with fluorescein isothiocyanate (FITC), and prepared from fresh plasma using chromatography. Plasminogen has been depleted by lysine affinity chromatography.
MOLECULAR WEIGHT:	330 kDa
CONCENTRATION:	8.43 mg/mL
STORAGE BUFFER:	0.02M Sodium Citrate; pH 7.4
STORAGE:	-70°C
PURITY:	>95% by SDS-PAGE analysis
SOURCE:	Rabbit Plasma
USAGE:	Thaw fibrinogen in a 37°C water bath. Aliquot as desired and freeze until needed. Keep at 25-37°C while using. Mix gently by pipet prior to use. Do not vortex.