

Product Specification Sheet

Biotin Labeled Human
Multimeric Vitronectin
Product Number: GP32
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Prepared from fresh human plasma using urea as a denaturant, then biotin labeled at primary amines.
MOLECULAR WEIGHT:	75,000 Monomer
CONCENTRATION:	2.5 mg/mL
SIZE:	0.1 mg
VOLUME:	0.04 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Human plasma
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; pH 7.4

For Research Purposes Only.