

**Product Specification
Sheet**

Porcine Vitronectin
Product Number: GP33
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Prepared from fresh porcine (pig) plasma using urea as a denaturant.
MOLECULAR WEIGHT:	75,000 Monomer
CONCENTRATION:	5 mg/mL
SIZE:	0.1 mg
VOLUME:	0.02 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Porcine plasma
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; pH 7.4

For Research Purposes Only.