

Product Specification Sheet

Anti-Human GST P1-1 Antisera
Product Number: GS72
Aliquot: 100 μ L
Lot Number: 031596DM
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum and is known not to cross-react with the mu or alpha class GST's.
MOLECULAR WEIGHT:	24 kDa
CONCENTRATION:	Not Applicable
STORAGE BUFFER:	None; Contains 0.05% sodium azide.
STORAGE:	-70°C. Aliquot once thawed and stored at -70°C. Avoid repeated freeze thaw cycles.
IMMUNOGEN:	Recombinant Human GST P1-1 expressed in <i>E. coli</i> .
PURIFICATION:	Not Purified
PURITY:	Not Purified
APPLICATIONS:	Western blot: 1:1,000 Is reactive with ≥ 0.2 μ g of recombinant human GST P1-1 (product # GS 70) in conjunction with an HRP conjugated donkey anti-rabbit secondary antibody.
SPECIES SPECIFICITY:	Human and Rat
REFERENCE:	Holmes, W., et al. <i>Br. J. Cancer</i> , 62 : 209-212, (1990)