

Product Specification Sheet

Anti-Human CYP450 3A4 IgG
Product Number: PA32
Aliquot: 100 μ L
Lot Number: pa32.091001
Storage: -70°C

- DESCRIPTION:** This polyclonal antibody is produced in rabbit and is provided as antiserum.
- MOLECULAR WEIGHT:** ~56 kDa
- STORAGE BUFFER:** TBS
- STORAGE:** -70°C. AVOID REPEATED FREEZE THAWS
- IMMUNOGEN:** Recombinant Human Cytochrome P450 3A4
- PURIFICATION:** None
- CLASS:** Antiserum
- APPLICATIONS:** Western blot: 1:2000 yields a positive result with 0.2 μ g of recombinant cytochrome P450 3A4 using HRP conjugated anti-rabbit secondary antibody.
- SPECIES SPECIFICITY:** Human.
- REFERENCES:** Gillam, et. al, Arch. Biochem & Biophys. 305: 123-131, 1993.

