

**Product Specification Sheet**

Anti-Human tPA Biotin Labeled  
IgG  
Product Number: PA33  
Aliquot: 1.0 mg  
Lot Number: pa33.1106  
Storage: -70°C

|                             |  |
|-----------------------------|--|
| <b>DESCRIPTION:</b>         | This rabbit polyclonal antibody is specific for human tissue plasminogen activator (tPA). The IgG preparation is provided with a biotin label. |
| <b>MOLECULAR WEIGHT:</b>    | 160 kDa  |
| <b>CONCENTRATION:</b>       | 7.3 mg/mL determined spectrophotometrically.   |
| <b>STORAGE BUFFER:</b>      | 0.5 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6  |
| <b>STORAGE:</b>             | -70°C. <b>AVOID REPEATED FREEZE THAWS.</b>   |
| <b>IMMUNOGEN:</b>           | Recombinant Human tPA  |
| <b>CLASS:</b>               | IgG  |
| <b>APPLICATIONS:</b>        | This antibody is suitable for use in Western blotting and ELISA application.<br>Western: 1:300 dilution<br>ELISA: 1:1000 dilution              |
| <b>SPECIES SPECIFICITY:</b> | Human  |
| <b>REFERENCES:</b>          | Spence, MJ <i>et al.</i> <i>Cytokine</i> <b>18:1</b> pp 26-34 (April 7, 2002)  |