

Product Specification Sheet

Sheep Anti-Mouse & Rat tPA
Antiserum
Product Number: PA38
Aliquot: 1.0 mL
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Anti-tPA (Mouse & Rat) Sheep Polyclonal Antiserum
MOLECULAR WEIGHT:	N/A
CONCENTRATION:	Whole antiserum
SIZE:	1.0 mL
VOLUME:	1.0 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Sheep serum
BUFFER:	N/A

For Research Purposes Only.