

Product Specification Sheet

Human sc-tPA, ATA-FPR-CMK
Inactivated
Product Number: PA42
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Prepared from human tPA by incubation with an excess of ATA-Phe-Pro-Arg-CMK. After labeling at the active site, remaining free peptide is removed by dialysis.
MOLECULAR WEIGHT:	70,000
CONCENTRATION:	1.0 mg/mL
SIZE:	0.1 mg
VOLUME:	0.1 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Chemical
BUFFER:	0.4 M HEPES; 0.1 M NaCl; pH 7.4

For Research Purposes Only.