

Anti-Mouse tPA MoAB
Product Number: PA62
Aliquot: 1.0 mg
Lot Number: pa62.1003A
Storage: -70°C

DESCRIPTION:	This is a capture monoclonal antibody produced in a tPA knockout mouse.
MOLECULAR WEIGHT:	70 kDa
CLONE:	H27B6
CONCENTRATION:	1.1 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.4
STORAGE:	-70°C
IMMUNOGEN:	Recombinant Mouse tPA
PURIFICATION:	Purified by Protein A Sepharose.
CLASS:	IgG
PURITY:	≥ 98% assessed by visual inspection of a Coomassie® blue-stained SDS-PAGE gel.
APPLICATIONS:	Suitable for use in ELISA, Western blot and immunohistochemistry. Recommended Western concentration at 1-5 µg/mL.
SPECIES SPECIFICITY:	Mouse and rat.