

Product Specification Sheet

Anti-Human uPA Affinity
Purified
Product Number: PA69
Aliquot: 0.1 mg
Lot Number pa69.706
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is specific for human uPA. The antibody reacts with both free and complex forms of uPA.
MOLECULAR WEIGHT:	54 kDa High Molecular Weight form 31 kDa Low Molecular Weight form
CONCENTRATION:	0.75 mg/mL determined Spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1m M EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant Human uPA
CLASS:	IgG
PURIFICATION:	Affinity purified by immobilized uPA.
APPLICATIONS:	Suitable for use in ELISA and Western blot.
SPECIES SPECIFICITY:	Human