

**Product Specification Sheet**

Mouse tPA, Biotin Labeled  
Product Number: PA72  
Aliquot: 0.05 mg  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Recombinant, produced in insect cells and biotin labeled.
<b>MOLECULAR WEIGHT:</b>	70,000
<b>CONCENTRATION:</b>	2.0 mg/mL
<b>SIZE:</b>	0.05 mg
<b>VOLUME:</b>	0.025 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Insect cell culture
<b>BUFFER:</b>	0.4M HEPES; 0.1M NaCl; pH 7.4

For Research Purposes Only.