

Product Specification Sheet

Recombinant Active Rabbit tPA
Product Number: PA77
Aliquot: 0.1 mg
Lot Number: pa77.210
Storage: -70°C

DESCRIPTION:	Recombinant Rabbit tPA produced in insect cells as a stable clone.
MOLECULAR WEIGHT:	70,000
CONCENTRATION:	1.3 mg/mL
SIZE:	0.1 mg
VOLUME:	0.077 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Insect cell culture; stable clone
ACTIVITY:	>95%
SPECIFIC ACTIVITY:	2,105,900 IU/mg (WHO standards for human tPA; NIBSC 98/714)
STORAGE BUFFER:	0.5M HEPES; 0.5M NaCl; pH 7.4

For Research Purposes Only.