

Product Specification Sheet

Recombinant Active Rat tPA
Product Number: PA78
Aliquot: 0.1 mg
Lot Number: pa78.110
Storage: -70°C

DESCRIPTION:	Recombinant Rat tPA produced in insect cells. >85% single chain and > 95% active
SPECIFIC ACTIVITY	1,060,000 IU/mg
MOLECULAR WEIGHT:	70,000 Da
CONCENTRATION:	4.88 mg/mL
SIZE:	0.1 mg
VOLUME:	0.02 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Insect cell culture
BUFFER:	0.4 M HEPES; 0.1 M NaCl; pH 7.4

For Research Purposes Only.