



### Product Specification Sheet

Description: Human UPA-FITC Labeled  
Product Number: PA81  
Aliquot: 0.1 mg  
Lot Number: pa81.617  
Storage: -70°C

**DESCRIPTION:** Human HMW Urokinase, Fluorescein Labeled

**MOLECULAR WEIGHT:** 54 kDa

**CONCENTRATION:** 1.85 mg/mL; 0.28 mole FITC / mole uPA

**STORAGE BUFFER:** 0.05 M Sodium Acetate; 0.1 M NaCl; 1 mM EDTA; pH 5.0

**STORAGE:** Store at -70°C and protect from light.

**SOURCE:** Human Urine

**PURITY:** >95% by SDS PAGE analysis