

Product Specification Sheet

Description: Amino Terminal
Fragment of uPA
Product Number: PA84
Aliquot: 0.023ml
Lot Number: pa84.711
Storage: -70°C

DESCRIPTION:	Amino terminal fragment of human urokinase purified from auto-catalytic digestion of recombinant human urokinase produced in insect cells. This fragment contains the Kringle and Epidermal Growth Factor domains.
MOLECULAR WEIGHT:	18,500
CONCENTRATION:	4.3 mg/ml
STORAGE BUFFER:	0.1M Tris-HCl; 0.15M NaCl; pH 7.4
STORAGE:	-70°C
SOURCE:	Insect Cell Culture
PURITY:	>95% by SDS-PAGE Analysis
EXTINCTION COEFFICIENT:	1.6