

**Product Specification
Sheet**

Recombinant Human
Cytochrome *b*₅
Product Number: PH91
Aliquot: 100 µg
Lot Number ph91.130510
Storage: -70°C

DESCRIPTION:	Active recombinant human cytochrome <i>b</i> ₅ protein. Histidine tagged.
CONCENTRATION:	1.77 mg/mL total protein using the BCA protein assay with BSA as a standard.
STORAGE BUFFER:	10 mM Na/K phosphate, pH 7.25, 0.1 mM EDTA, 0.1 mM DTT, 0.01% CHAPS
STORAGE:	Store at -70°C. Avoid Multiple Freeze-Thaw Cycles.
PURITY:	≥75% by visual inspection of 5 µg in a Coomassie® stained SDS gel.
SOURCE:	<i>E. coli</i> expressed histidine-tagged recombinant human protein.
REFERENCE:	Holmans <i>et al.</i> (1994) <i>Arch. Biochem. Biophys.</i> 312 : 554-565.