

**Product Specification Sheet**

Human PAI-1 Latent Wild Type  
Product Number: PI02  
Aliquot: 0.5 mg  
Lot Number: pi02.1109  
Storage: -70°C

**DESCRIPTION:** This preparation is <0.4% active as determined by titration with HMW tc-urokinase and SDS PAGE. Latent PAI-1 is a useful control for experiments when used in conjunction with the active form above.

**MOLECULAR WEIGHT:** 43 kDa

**CONCENTRATION:** 1.9 mg/mL


**STORAGE BUFFER:** 0.05 M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6

**STORAGE:** -70°C. Avoid repeated Freeze-Thaw cycles.

**SOURCE:** E. coli

**FORM:** Frozen liquid

**PURITY:** ≥95% by SDS-PAGE gel.



10% SDS-PAGE - GelCode Blue Stain  
1 PAI-L-697 (5 ug) Reduced  
2 PAI-L-697 (5 ug)+uPA-HTC (8 ug) Reduced  
3 Prestained Standard

For in vitro research use only.