

**Product Specification Sheet**

Human PAI-1 (Stable Mutant)  
Product Number: PI05  
Aliquot: 0.5 mg  
Lot Number pi05.320  
Storage: -70°C

<b>DESCRIPTION:</b>	Human Plasminogen Activator Inhibitor 1 (Stable Mutant)
<b>CONCENTRATION:</b>	2.2 mg/mL
<b>MOLECULAR WEIGHT:</b>	43,000
<b>STORAGE BUFFER:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>STORAGE:</b>	Store at -70°C. Avoid multiple freeze/thaw cycles.
<b>VOLUME:</b>	0.5 mL
<b>SPECIFIC ACTIVITY:</b>	>95% active by uPA titration.
<b>SOURCE:</b>	E. coli