

Product Specification Sheet

Human PAI-1, His tagged
Product Number: PI06
Aliquot: 0.5 mg
Lot Number pi06.1010
Storage: -70°C

DESCRIPTION:	Human Plasminogen Activator Inhibitor-1, Histidine tagged (stable mutant).
QUANTITY:	0.5 mg
MOLECULAR WEIGHT:	43,000
STORAGE:	This product is to be stored at -70°C and is stable for one year from the date of arrival.
CONCENTRATION:	1.9 mg/mL
VOLUME:	0.263 mL
BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6.