

**Product Specification Sheet**

Rb anti-Human PAI-1 antibody  
Product Number: PI34  
Aliquot: 0.1 mg  
Lot Number: pi34.810  
Storage: -70°C

<b>DESCRIPTION:</b>	Rabbit polyclonal antibody against human PAI-1, Affinity purified by immobilized human PAI-1.
<b>MOLECULAR WEIGHT:</b>	43,000
<b>CONCENTRATION:</b>	4.8 mg/mL
<b>STORAGE BUFFER:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Recombinant full length wt- human PAI-1
<b>PURIFICATION:</b>	Affinity purified
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	Suitable for ELISA, IHC, and WB