

**Product Specification Sheet**

Sheep Anti-Human PAI-1  
Antiserum  
Product Number: PI57  
Aliquot: 1.0 mL  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Anti-PAI-1 (human), sheep polyclonal antiserum.
<b>MOLECULAR WEIGHT:</b>	N/A
<b>CONCENTRATION:</b>	Whole antiserum
<b>SIZE:</b>	1.0 mL
<b>VOLUME:</b>	1.0 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Sheep serum
<b>BUFFER:</b>	N/A

For Research Purposes Only.