

Product Specification Sheet

Anti-2,8,17 Citrullinated Histone
H3

Product Number: NF24

Aliquot: 100 µg

Lot Number: nf24.220306

Storage: -20°C

DESCRIPTION:	This mouse monoclonal antibody to human 2,8,17-Citrullinated Histone H3 is protein G affinity purified.
CONCENTRATION:	1 mg/mL
STORAGE BUFFER:	Phosphate Buffered Saline pH 7.4, 50% Glycerol, 0.02% Proclin 150
STORAGE:	-20°C. Avoid repeated freeze thaw cycles
IMMUNOGEN:	Synthetic peptide coupled to KLH.
PURIFICATION:	Protein G Affinity Purification
APPLICATIONS:	ELISA, Western Blotting
SPECIES SPECIFICITY:	Human