

Product Specification Sheet

Anti-Troponin T
Product Number: NF35
Aliquot: 100 µg
Lot Number: nf35.220308
Storage: -20°C

DESCRIPTION: This mouse monoclonal antibody to human Troponin T is protein G affinity purified.

CONCENTRATION: 1 mg/mL

STORAGE BUFFER: Phosphate Buffered Saline pH 7.4, 50% Glycerol, 0.02% Proclin 150

STORAGE: -20°C. Avoid repeated freeze thaw cycles

IMMUNOGEN: Synthetic Peptide Coupled to KLH

PURIFICATION: Protein G Affinity Purified

APPLICATIONS: ELISA, Western Blotting

SPECIES SPECIFICITY: Human