

Product Specification Sheet

Anti-VEGF
Product Number: NF49
Aliquot: 100 µg
Lot Number: nf49.220311
Storage: -20°C

DESCRIPTION:	This mouse monoclonal antibody to human VEGF is protein G affinity purified.
CONCENTRATION:	1 mg/mL
STORAGE BUFFER:	Phosphate Buffered Saline pH 7.4, 50% Glycerol, 0.02% Proclin 150
STORAGE:	-20°C. Avoid repeated freeze thaw cycles
IMMUNOGEN:	Recombinant Human VEGF expressed in <i>E.coli</i> .
PURIFICATION:	Protein G Affinity Purified
APPLICATIONS:	ELISA, Western Blotting
SPECIES SPECIFICITY:	Human