

Product Specification Sheet

Anti-NT ProBNP
Product Number: NF56
Aliquot: 100 µg
Lot Number: nf56.220310
Storage: -20°C

DESCRIPTION: This mouse monoclonal antibody to human NT-ProBNP is protein G affinity purified.

CONCENTRATION: 1 mg/mL

STORAGE BUFFER: Phosphate Buffered Saline pH 7.4, 50% Glycerol, 0.02% Proclin 150

STORAGE: -20°C. Avoid repeated freeze thaw cycles

IMMUNOGEN: Synthetic Peptide Coupled to KLH

PURIFICATION: Protein G Affinity Purification

APPLICATIONS: ELISA, Western Blotting

SPECIES SPECIFICITY: Human