

**Product Specification Sheet**

Anti-Mouse PAI-1 (MoAB) -  
Capture  
Product Number: PI67  
Aliquot: 1.0 mg  
Lot Number pi67.1211  
Storage: -70°C

<b>DESCRIPTION:</b>	This capture monoclonal antibody is produced in a PAI-1 knockout mouse.
<b>MOLECULAR WEIGHT:</b>	43 kDa non-glycosylated 50 kDa glycolysated
<b>CLONE:</b>	H34G6
<b>CONCENTRATION:</b>	7.6 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 Sodium Phosphate; 0.15 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Mouse PAI-1
<b>PURIFICATION:</b>	Affinity purified by immobilized PAI-1
<b>PURITY:</b>	≥ 90% assessed by visual inspection of a Coomassie® blue-stained SDS-PAGE gel.
<b>CLASS</b>	IgG
<b>APPLICATIONS:</b>	ELISA
<b>SPECIES SPECIFICITY:</b>	Mouse