

Product Specification Sheet

Anti-Human Protein Kinase C
gamma (PKC γ)
Product Number: PK29
Aliquot: 100 μ L
Lot Number: pk29.typical lot
Storage: -20° C

IMMUNOGLOBULIN:	Polyclonal Rabbit IgG
IMMUNOGEN:	PKC γ Peptide Sequence: ³⁰² F-E-A-C-N-Y-P-L-E-L-Y-E-V-S-R-Y-G-P-S-S-S ³²²
SPECIFICITY:	Prepared against synthetic peptide corresponding to distinct C terminal region. It is based on known homologies and is expected to react with PKC orthologs in most mammalian species.
STORAGE:	Store at -20° C. Avoid repeating freeze/thaw cycles, as it may result in loss of activity.
PURITY:	Serum
WORKING DILUTION:	Suggested working dilution for Western blots is 1:1000. For immunohistochemistry, a dilution of 1:250 is recommended. This product may also be used in ELISA applications.
REFERENCES:	Wetsel, W.C., et al., <i>J. Cell Biol.</i> 117: 121-133 (1992)