

Anti-Protein Kinase C ϵ
Product Number: PK35
Aliquot: 100 μ L
Lot Number: 102302CS
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C ϵ (PKC ϵ).
MOLECULAR WEIGHT:	83.7 kDa Calculated
CONCENTRATION:	Not Determined
ISOTYPE:	NA
STORAGE BUFFER:	Phosphate Buffered Saline; 0.05% Sodium Azide
STORAGE:	-70°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
IMMUNOGEN:	KLH conjugated synthetic peptide corresponding to amino acids 721- 737 of PKC ϵ . Sequence: ⁷²¹ N-Q-E-E-F-K-G-F-S-Y-F-G-E-D-L-M-P ⁷³⁷
PURIFICATION:	NA
APPLICATIONS:	Western Blotting: 1:1,000 Dilution
SPECIES SPECIFICITY:	Human, mouse and rat among other mammalian species.
REFERENCE:	JBC (1993) 268 :20110-20115