

Product Specification Sheet

Description: Protein Kinase C Eta
Product Number: PK46
Aliquot: 1 µg
Lot Number: 104693JC
Storage: -70°C

DESCRIPTION:	Human recombinant kinase C η expressed in baculoviral insect expression system. The protein is suspended in SDS/mercaptoethanol sample buffer with a bromophenol blue dye marker.
RECOMMENDED DILUTION:	Load 2-5 µL of prepared standard directly from vial. No dilution is necessary.
CONCENTRATION:	20 µg/mL.
ENZYME ACTIVITY:	Not enzymatically active.
STORAGE:	Stable indefinitely at -70°C; for at least one year at -20°C.
APPLICATIONS:	Intended for use as a standard for Western blots using enhanced chemiluminescence. For blots using less sensitive substrates, such as alkaline phosphatase, larger amounts of protein should be loaded onto the gel.
PACKAGING:	Packaged as 1 µg per vial of protein prepared in sample buffer.