

**Product Specification Sheet**

Anti-Protein Kinase C theta ( $\theta$ )  
Product Number: PK55  
Aliquot: 100  $\mu$ L  
Lot Number 020296YY  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C $\theta$ (PKC $\theta$ ).
<b>CONCENTRATION:</b>	Not Determined
<b>PRESENTATION:</b>	Antiserum with 0.05% Sodium Azide
<b>STORAGE:</b>	-70°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
<b>IMMUNOGEN:</b>	KLH conjugated synthetic peptide corresponding to amino acids found on the C-terminus of PKC $\theta$ . <b>Sequence:</b> H-S-N-F-S-F-I-N-P-G-M-E-T-L-I-C-S
<b>PURIFICATION:</b>	Not Purified
<b>APPLICATIONS:</b>	Western Blotting: 1:250dilution using alkaline phosphatase detection technology.
<b>SPECIES SPECIFICITY:</b>	Human, mouse and pig among other mammalian species.