

Product Specification Sheet

Human Glutamine Plasminogen
Product Number: PL20
Aliquot: 1.0 mg
Lot Number: p120.716
Storage: -70°C

DESCRIPTION:	This is human glutamic acid plasminogen. It is prepared from plasma from lysine sepharose chromatography.
MOLECULAR WEIGHT:	92 kDa
CONCENTRATION:	9.8 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.1 M HEPES; 0.1 M NaCl; pH 7.4
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
PURIFICATION:	Prepared by immobilized lysine chromatography.
ACTIVITY:	No plasmin activity is detected with chromogenic substrate S-2251. >98% conversion to plasmin is observed upon activation with uPA.