

**Product Specification Sheet**

Rabbit Plasmin  
Product Number: PL38  
Aliquot: 1.0 mg  
Lot Number: pl38.317  
Storage: -70°C

<b>DESCRIPTION:</b>	Rabbit plasmin prepared from plasminogen by activation with immobilized human uPA
<b>MOLECULAR WEIGHT:</b>	85 kDa
<b>CONCENTRATION:</b>	8.3 mg/mL
<b>STORAGE BUFFER:</b>	0.1 M HEPES; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAW CYCLES.</b>
<b>SOURCE:</b>	Rabbit plasma
<b>PURITY:</b>	≥ 95 % by SDS PAGE
<b>NOTES:</b>	Plasmin will undergo rapid autolysis in the absence of protease inhibitors such as benzamidine, and should be used quickly after thawing before significant activity losses occur.