

Product Specification Sheet

Anti-Mouse Antiplasmin IgG
Fraction
Product Number: PL49
Aliquot: 1.0 mg
Lot Number: p149.506A
Storage: -70°C

DESCRIPTION: This rabbit polyclonal antibody is specific for mouse antiplasmin and is provided as an IgG preparation.

MOLECULAR WEIGHT: 67 kDa for mature antiplasmin**

CONCENTRATION: 2.1 mg/mL determined spectrophotometrically.

STORAGE BUFFER: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6

STORAGE: -70°C. **AVOID REPEATED FREEZE THAWS.**

IMMUNOGEN: Mouse Antiplasmin

ISOTYPE: IgG

PURIFICATION: Purified by Protein A Sepharose.

APPLICATIONS: Western blotting

SPECIES SPECIFICITY: Mouse

****REFERENCE:** *Hanss, Michel ML a et al., Blood Coagulation & Fibrinolysis. 14(1):107-111, January 2003.*