

Product Specification Sheet

Anti-Mouse Antiplasmin MoAB
IgG (Inhibitory)
Product Number: PL57
Aliquot: 1.0 mg
Lot Number: p157.804A
Storage: -70°C

DESCRIPTION: This mouse monoclonal antibody produced in an antiplasmin knockout mouse. The antibody is provided as an IgG prep and is inhibitory.

MOLECULAR WEIGHT: 67 kDa for mature antiplasmin**

CONCENTRATION: 1.9 mg/mL determined spectrophotometrically.

CLONE: MAP27C9

STORAGE BUFFER: 0.5 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.8

STORAGE: -70°C. **AVOID REPEATED FREEZE THAWS.**

IMMUNOGEN: Mouse Antiplasmin

ISOTYPE: IgG

PURIFICATION: Purified by Protein A Sepharose.

APPLICATIONS: Western blotting

SPECIES SPECIFICITY: Mouse

****REFERENCE:** *Hanss, Michel ML a et al.*, Blood Coagulation & Fibrinolysis. 14(1):107-111, January 2003.