

Product Specification Sheet

Anti-Murine Plasminogen/Plasmin
Anti-Serum
Product Number: PL60
Aliquot: 1.0 mg
Lot Number: pl60.905A
Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is generated against mouse plasminogen and will cross react with plasmin. The antibody is provided as anti-serum.
CONCENTRATION:	Not Applicable
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Mouse Plasminogen
PURIFICATION:	Not Purified
APPLICATIONS:	Western blotting and ELISA
SPECIES SPECIFICITY:	Mouse