



Product Specification Sheet

Rabbit Anti-Murine Plasminogen
(IgG Fraction)
Product Number: PL62
Aliquot: 1.0 mg
Lot Number: pl62.506
Storage: -70°C

DESCRIPTION:	Rabbit Anti-Murine Plasminogen (IgG Fraction)
CONCENTRATION:	11.2 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Mouse Plasminogen
ISOTYPE:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
APPLICATIONS:	Western blot: 1:1000 dilution ELISA: 1:10,000-1:100,000 dilution
SPECIES SPECIFICITY:	Mouse.