

**Product Specification Sheet**

Anti-Mouse Plasminogen Affinity  
Purified IgG  
Product Number: PL64  
Aliquot: 0.1 mg  
Lot Number: pl64.105A  
Storage: -70°C

<b>DESCRIPTION:</b>	This rabbit polyclonal antibody is specific for mouse plasminogen and is useful in cancer and fibrinolytic studies. This is a high titer preparation.
<b>MOLECULAR WEIGHT:</b>	92 kDa
<b>CONCENTRATION:</b>	2.9 mg/mL determined spectrophotometrically.
<b>PURITY:</b>	> 98% pure by visual inspection on a Coomassie™ Blue stained SDS-PAGE.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Purified Mouse Plasminogen
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Affinity purified by immobilized plasmin.
<b>APPLICATIONS:</b>	This preparation is suitable for use in Western blotting and ELISA.
<b>SPECIES SPECIFICITY:</b>	Mouse